EN - Examination gloves

Synthetic polymer (PVC-Nitrile) * Latex Free * Powder Free * Non Sterile * Ambidextrous
* Warning: Gloves should be kept in their original packaging * Store in a cool dry place * Keep away from sources of ozone or ignition

Personal Protective Equipment of Complex Design Category III, in compliance with Regulation (EU) 2016/425, type tested to EN 420:2003+A1:2009, EN ISO 374-1:2016+A1:2018, EN 374-2:2014, EN 16523-1:2015, EN 374-4:2013 and EN ISO 374-5:2016, CE 2797; Class I Medical Device, in compliance with MDR 2017/745, type tested to 455-1:2020, 455-2:2015, 455-3:2015, EN 455-4:2009 and ASTM D5250.

This product provides protection against Bacteria, Fungi and Virus. The gloves had been tested in accordance with ISO 16604:2004 to meet the requirements of EN ISO 374-5:2016 for resistance to penetration by blood-borne pathogens-test method using Phi-X174 bacteriophage. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.
Protection against bacteria and fungi - Pass
Protection against viruses - Pass

Notified Body 2797:
BSI Group The Netherlands BV, John M. Keynesplein 9, 1066 EP Amsterdam, NL

User information:

User: These gloves are designed for single use only and should be replaced after use. A liquid proof design allows for protection against certain chemicals. Following suitable risk assessment, typical applications may include: medical examination, healthcare, domestic and household, product examination, cleaning, food processing and industrial applications.
Precaution For Use: Before use, inspect for any defects or imperfections. If in doubt, discard the gloves and select a new pair.

If gloves are to be used with chemicals:

* Consult the subsequent user instruction to ensure these gloves are suitable for the intended purpose.
* Keep all chemicals away from direct skin contact, even if thought to be harmless. If contact occurs, wash the affected area immediately and seek assistance, contacting a suitable, qualified professional.
* Gloves that have been in contact with chemicals should be removed and disposed of as quickly as possible following use.
* Ensure chemicals or residuals cannot enter via the cuff and come in contact with skin.

* These gloves should not be used in applications which require mechanical and/or thermal protection.

* Discontinue use immediately if signs of tearing, swelling or degradation appear.

Ingredients/Hazardous ingredients: Some gloves may contain ingredients which are known to be a possible cause of allergies to sensitive individuals, including irritation and/or allergic reactions. If this occurs, seek appropriate medical advice immediately.

Storage: Keep out of direct sunlight, storing in a cool, dry, well ventilated and clean location. Keep away from sources of ozone or ignition. Gloves should be kept in their original packaging.

Cleaning: These gloves are not designed to be laundered or cleaned and are for single use only.

Disposal: Used gloves may be contaminated with infectious or other hazardous materials. Dispose of according to local authority regulations. Landfill or incinerate under controlled conditions.

Composition: This product is made from Synthetic polymer (PVC-Nitrile) and does not contain natural rubber.

Limitation/ Caution/ Warning

* The information enclosed is for products in a new condition and does not reflect the actual duration of protection in the workplace nor does it reflect the potential differentiation between 'mixtures' and 'pure chemicals'.

* Before use, it is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. Before usage, inspect the gloves for any defect or imperfections.

* The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only.

* When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.

* The glove length is 240 mm and is designed for protection till the cuff of hand and is an examination glove.

* The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

EN ISO 374-1:2016+A1:2018 permeation levels are based on breakthrough times as follows:

Performance Level	1	2	3	4	5	6
Measured breakthrough time (mins)	>10	>30	>60	>120	>240	>480

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