## EN - Examination gloves

Latex Free \* Powder Free \* Non Sterile \* Ambidextrous \* Warning: Gloves should be kept in their original packaging \* Store in a cool dry place \* Keep away from sources of ozone or ignition

Personal Protective Equipment Category III, in compliance with Regulation (EU) 2016/425, type tested to EN420:2003+A1:2009, EN ISO 374-1:2016+A1:2018, EN 374-2:2014, EN 16523-1:2015, EN 374-4:2013 and EN ISO 374-5: 2016, CE 2777; Class I Medical Device, in compliance with MDR 2017/745, type tested to EN 455-1:2020, EN 455-2:2015, EN 455-3:2015, EN 455-4:2009.

This product provides protection against Bacteria, Fungi and Virus. The gloves had been tested in accordance with ISO 16604:2014 to meet the requirements of EN ISO 374-5:2016 for resistance to penetration by blood-borne pathogens-test method using Phi-X174 bacteriophage. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.  
Protection against bacteria and fungi - Pass  
Protection against viruses - Pass

Notified Body responsible for EU Type Examination (Module B) and Module C2 on-going Conformity: 2777: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee D15 YN2P, Ireland

### User information:

**Use :** These gloves are designed for single use only and should be replaced after use. A liquid proof design allows for protection against certain chemicals. Following suitable risk assessment, typical applications may include: medical examination, healthcare, domestic and household, product examination, cleaning, food processing and industrial applications.

**Precaution For Use:** Before use, inspect for any defects or imperfections. If in doubt, discard the gloves and select a new pair.

If gloves are to be used with chemicals:

\* Consult the subsequent user instruction to ensure these gloves are suitable for the intended purpose.

\* Keep all chemicals away from direct skin contact, even if thought to be harmless. If contact occurs, wash the affected area immediately and seek assistance, contacting a suitable, qualified professional.

\* Gloves that have been in contact with chemicals should be removed and disposed of as quickly as possible following use.

\* Ensure chemicals or residuals cannot enter via the cuff and come in contact with skin.

\* These gloves should not be used in applications which require mechanical and/or thermal protection.

\* Discard use immediately if signs of tearing, swelling or degradation appear.

**Ingredients/Hazardous ingredients:** Some gloves may contain ingredients which are known to be a possible cause of allergies to sensitive individuals, including irritation and/or allergic reactions. If this occurs, seek appropriate medical advice immediately. **The product contains dithiocarbamates - not to be used in the case of hypersensitivity to these substances.**

**Storage:** Keep out of direct sunlight, storing in a cool, dry, well ventilated and clean location. Keep away from sources of ozone or ignition. Gloves should be kept in their original packaging.

**Cleaning:** These gloves are not designed to be laundered or cleaned and are for single use only.

**Disposal:** Used gloves may be contaminated with infectious or other hazardous materials. Dispose of according to local authority regulations. Landfill or incinerate under controlled conditions.

**Composition:** This product is made from 100% synthetic nitrile polymer and does not contain natural rubber

### Limitation/ caution/ warning

\* The information enclosed is for products in a new condition and does not reflect the actual duration of protection in the workplace nor does it reflect the potential differentiation between 'mixtures' and 'pure chemicals'.

\* Before use, it is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. Before usage, inspect the gloves for any defect or imperfections.

\* The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture.

\* When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.

\* The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

Gloves for Special Applications (EN 420:2003+A1:2009, Clause 5.1.3). These gloves are designed to protect the hand surface from chemical splashes. Therefore, the length of the gloves is below EN requirements of total minimum length, and deems as 'Fit for Special Purpose'.

EN ISO 374-1:2016+A1:2018 permeation levels are based on breakthrough times as follows:

Performance Level	1	2	3	4	5	6
Measured breakthrough time (mins)	>10	>30	>60	>120	>240	>480

Test chemical	EN ISO 374-1:2016+A1:2018 Permeation Level	EN 374-4:2013 Degradation (mean value)
K Sodium Hydroxide 40%	6	- 9,5%
P Hydrogen peroxide 30%	6	+ 44,0%
T Formaldehyde 37%	4	+ 51,0%

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