




**MEDICA**  
EUROPE


# Medica Brush®

A disposable scrub brush with sponge



## Medica Brush® is designed for cleaning:

- nails with the nail cleaner
- nails and knuckles with the sponge
- hands and forearms with the sponge



## A cross cut shows:

- the soft long hairs for mild cleaning of the knuckles
- the short hairs for perfect cleaning under the nails



## Squeeze and release:

- squeezing the brush ensures the sponge to absorb more water for optimal cleaning of hands and forearms



Manufacturing of the Medica Brush® is in The Netherlands and the process is fully automatic which ensures a constant quality in assembling brush and sponge.